Magnetic couplings, stirrer drives for chemical Reactors bmd and cyclone: Absolutely tight, safe and dependable

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 **büchiglasuster**[®]
switzerland Pilot Plant and Reactor Systems

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Peter-Henlein-Str. 20
D-85540 Haar b. München
Telefon 089/45 60 06 70
Telefax 089/45 60 06 80
info@c3-analysentechnik.de
www.c3-analysentechnik.de

BÜCHI – THE WAY TO GET RESULTS!

Magnetic couplings

Magnetkupplungen


The main difference of a magnetic stirrer drive and a conventional stirrer drive such as rubber gland, lip seal and mechanical seal is the fact, that the internal stirrer shaft is driven without any contact to the static seal. The external magnet and the internal magnet are virtually connected through a strong magnetic field. The external magnet is located outside the pressure chamber. The internal magnet, which is inside the pressure chamber, is connected to the stirrer shaft and rotates simultaneously with the motor driven external magnet. Buchi offers low to high torque magnetic couplings and stirrer drives, which fulfill high performance mixing also in high viscous media.

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Functional principle of a bmd magnetic coupling

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bmd 800

bmd 1200

For technical Details of standard bmd 075 to bmd 5400 please refer to the table on reverse side

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Magnetic drives with 75 / 300 Ncm torque (cyclone 075 – cyclone 300)

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Büchi cyclone 075 as well as the cyclone 300 are complete stirrer drives. They contain a magnetic coupling bmd 075 / bmd 300 with an electric motor. They can be installed on the reactor cover plate with an M36x1.5 thread.

The cyclone 300 controller enables stepless stirrer speed adjustment of the servo motor. Stirrer speed and torque can be transferred via a RS 232 interface to bds measuring / display unit and/or to the bls software. The stirrer speed can be set electrically by a 4...20 mA Signal.

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