

**NEW**

**Cooling System for THINKY MIXER**

# ENS-10CE

Supported Models: THINKY MIXER ARE-250CE / Solder Paste Mixer SR-500CE

**The powerful DC fan can forcibly remove heat accumulated in the mixer system.**

Using ENS-10CE can lower operational temperature by 10 °C

- THINKY Cooling System ENS-10CE is effective in decreasing temperature by 10 °C in the mixer.
- THINKY Cooling System ENS-10CE is effective in suppressing temperature raises in the mixer caused by generation of heat by materials being mixed.
- Interval time for cool-down can be shortened and efficient repeat operation is possible.



Planetary Centrifugal Mixer  
**THINKY MIXER ARE-250CE**

Cooling System  
**ENS-10CE**

**On/Off button**  
 The ON and OFF button of a mixer power supply are also built into the ENS-10CE.

**Effect of ENS-10CE for repetitive operation**

**Example 1**  
**High Viscosity Silicone**  
 2000 rpm, 5 min  
 10 times repeatedly

**Example 2**  
**Medium Viscosity Polyimides**  
 2000 rpm, 15 min  
 15 times repeatedly

**Example 3**  
**Dispersing Phosphor into Resin**  
 400 rpm, 30 s on 1st step,  
 1500 rpm, a few min on 2nd step  
 10 times repeatedly



**SR-500CE**

\* The product appearance and specification may be changed without a preliminary announcement.

**ENS-10CE is able to suppress the raise in temperature inside of THINKY MIXER during repetitive operation.**

THINKY Cooling System ENS-10CE is a new optional cooling system that exhausts heat from THINKY MIXER. Temperatures in a mixer are 10 °C higher without the use of this cooling system.



**Test Report: Comparison table with and without using ENS-10CE**

Test material: Shin-Etsu Silicone, Silicone Fruide KF-96-5,000CS

**Results**

The effect was proved that temperatures in the mixer were 10 °C less and 11 °C less on the motor.

**DC fan OFF**



Temperature change of each part of the mixer

Power OFF	Average
Internal temperature	43.9 °C
Motor temperature	61.6 °C
DC exhaust temperature	34.1 °C
Room temperature	26.4 °C

**DC fan ON**

Temperature change of each part of the mixer

Power ON	Average
Internal temperature	33.2 °C
Motor temperature	50.4 °C
DC exhaust temperature	32.7 °C
Room temperature	25.7 °C



**Setting Image**



● ARE-250CE



● SR-500CE

**Product Specification**

Model	ENS-10CE	
Usage	Reduces the temperature in THINKY MIXER	
Supported models	ARE-250CE, SR-500CE	
Temperature reduction effects	Temperature in the mixer : Reduced approx. 10°C *1	
Power supply	Voltage	Single-phase AC230 V± 10 %, 50 Hz
	Power Consumption	Exhauster: about 40 VA mixer combination: Maximum 900 VA (operation)
Environment	10 to 35°C, 35 to 85% RH (thing without dew condensation)	
Dimension	H 145 mm × W 310 mm × D 320 mm	
Unit Weight	Approx. 7 kg	

\*1: Please contact us for higher specification.

The product appearance and specification may be changed without a preliminary announcement.

**THINKY CORPORATION**

Headquarters: 2-16-2 Sotokanda, Chiyoda-ku, Tokyo 101-0021  
Phone : +81-3-5207-2713 Fax : +81-3-5289-3281

THINKY USA Inc. : 23151 Verdugo Drive, Suite 112 Laguna Hills, CA92653, USA

THINKY CHINA: East building, HaiAn Kafunuo Mansion, Shennan road, Qianhai road, Nanshan district, Shenzhen

<https://www.thinkymixer.com/en-gl/>

**Sales agent**



C3 PROZESS- UND  
ANALYSENTECHNIK GmbH  
Peter-Henlein-Str. 20  
D-85540 Haar b. München  
Telefon 089/45 60 06 70  
Telefax 089/45 60 06 80  
info@c3-analysentechnik.de  
www.c3-analysentechnik.de