

The new Büchi pressure filters – filtration on a new level

- easy operation
- removable filter basket
- quick opening / closing system
- perfect to Büchi pressure reactors
- for every-day lab use



Facts and figures

- Made from stainless steel or Hastelloy
- Volume up to 20 litres
- Up to 6 bar pressure and vacuum operation
- Tri-Clamp connections
- Machinery directive/PED

The new pressure filter by Büchi – simple and superior

Accessories

- Various filter sizes and filter types available
- With or without heating jacket, adjusted to your needs



Small version up to a volume of 5 litres.



Large version up to a volume of 20 litres.



Clamp screws guarantee easy opening.



Simply lift the filter basket for removal of the filter cake and cleaning. A time-saving and cost efficient solution.



Visit www.buchiglas.com
for more information about:

- pressure reactors
- glass plants
- automation

 **büchiglasuster**
switzerland Pilot Plant and Reactor Systems

Büchi AG
Gschwaderstrasse 12
CH-8610 Uster/ Switzerland

Phone +41 (0) 44 905 51 11
Fax +41 (0) 44 905 51 22
www.buchiglas.com



C3 PROZESS- UND
ANALYSENTECHNIK

Peter-Henlein-Str. 20
D-85540 Haar b. München
Telefon 089/45 60 06 70
Telefax 089/45 60 06 80
info@c3-analysentechnik.de
www.c3-analysentechnik.de